










# VORTRAN Automatic Resuscitator (VAR)<sup>TM</sup>

Unique single patient, multiple-use disposable emergency resuscitator

<b>VAR-Plus Model PCM for patients between 10 - 40 kilograms</b>	<b>VAR Models RT and RC for patients 40 kilograms</b>		
	 <input type="checkbox"/> Model RT	 <input type="checkbox"/> RT w/Manometer	 <input type="checkbox"/> RT w/Extension
 <input type="checkbox"/> PCM w/Manometer	 <input type="checkbox"/> Model RC	 <input type="checkbox"/> RC w/Manometer	 <input type="checkbox"/> RC w/Extension

## USER'S GUIDE

Revised Rev: 4/14/2005

### TABLE OF CONTENTS

- Table of Contents ..... 1
- I. Functional and Operational Characteristics ..... 2
  - Figure 1 - VAR Component Description ..... 2
- II. Clinical Considerations ..... 3
  - Figure 2 - Airway Pressures - PIP & PEEP ..... 6
- III. Protocol: Setup Instructions - VAR ..... 8
- IV. Cautions and Warnings ..... 12
- V. VAR Competency ..... 13
- VI. Frequently Asked Questions ..... 14
- VII. Clinical Reference ..... 19
- VIII. Coding Information ..... 23
- IX. Troubleshooting ..... 24
- X. Ordering Information ..... 25
- XI. Quick Guide ..... 26


## VAR User's Guide

### I. Functional and Operational Characteristics

The gas-powered VAR provides constant flow, pressure cycled automatic ventilatory support for both breathing and non-breathing patients. The primary working mechanism of the VAR (refer to Figure 1) is the <sup>[1]</sup> modulator with <sup>[a]</sup> peak inspiratory pressure (PIP) and the <sup>[b]</sup> breathing rate adjustment dials, which includes an exhalation valve that opens at one pressure (PIP) and closes at another lower pressure (PEEP). The remaining components of the VAR consist of the <sup>[2]</sup> pressure manometer, <sup>[3]</sup> gas inlet for supply gas flow, <sup>[4]</sup> patient connection port, <sup>[5]</sup> redundant pressure pop-off valve, and <sup>[6]</sup> one-way valve for entraining additional air.

The pulmonary modulator provides the actual ventilatory support. The primary working mechanism of the pulmonary modulator is the piston. The piston is spring loaded, designed like a pressure pop-off valve except the spring force is adjustable (the <sup>[a]</sup> PIP Dial).

Figure 1  
VAR (Model RCM) Component Description



### OPERATIONAL CHARACTERISTICS

	<i>VAR-Plus</i> (Model PCM)	VAR (Model RT & RC)
Recommended body weight	10 kg and above	Greater than 40 kg (adults)
Ventilatory frequency	Auto-adjusting to lung capacity	8 to 20 BPM
Adjustable PIP range	10 to 45 cm-H <sub>2</sub> O	20 to 50 cm-H <sub>2</sub> O
PEEP	1/5 <sup>th</sup> of PIP (2 to 9 cm-H <sub>2</sub> O)	1/10 <sup>th</sup> of PIP (2 to 5 cm-H <sub>2</sub> O)
Inspiratory resistance	3 ± 1 cm-H <sub>2</sub> O/L/sec	SAME
Expiratory resistance	3 ± 1 cm-H <sub>2</sub> O/L/sec	SAME
Dead space	4 ± 3 mL	SAME
Operating environmental limits	-18 to 50°C	SAME
Storage environmental limits	-40 to 60°C	SAME
Patient connection	15 mm female, 22 mm male	SAME
Gas inlet	DISS connection	SAME
Oxygen concentration	50-85% when supplied w/100% O <sub>2</sub>	SAME

## VAR User's Guide

---

### II. Clinical Considerations

The VAR provides short term, pressure cycled, and constant flow ventilatory support for either breathing or non-breathing patients. This allows the patient to receive consistent and reliable ventilatory support. Because the VAR is pressure cycled, changes in the patient's lung compliance will cause a change in the patient's breathing rate. The VAR is positional sensitive. Final adjustments should be made with the VAR in its secured operating position. The VAR is pressure cycled on inhalation and exhalation (PIP and PEEP) which minimizes the possibility of gas trapping. During inhalation, exhalation will not start until PIP is reached. During exhalation, inhalation will not begin until pressure drops to PEEP. For the spontaneously breathing patient, the rate dial of the VAR is set so that the base line pressure is above the intrinsic PEEP allowing the patient to initiate inhalation by drawing the base line pressure down to the set PEEP. Because the VAR is a constant flow pressure cycled device, changes in patient compliance will result in changes in the respiratory rate (stiffer or smaller compliances produce faster rates). The advantage of this minimizes the danger of barotrauma. It should be emphasized that the VAR is to be used only by trained personnel who continuously monitor the patient. The VAR is not an ICU stand-alone ventilator with multiple monitoring features.

Setup and use of the VAR is simple (refer to Setup Instructions in Section III on page 7). Set desired flow ( $Q$ ), adjust PIP pressure dial to obtain desired inspiratory time ( $t_{insp}$ ) to attain tidal volume ( $TV = Q \times t_{insp}$  see Tidal Volume Table 1). The gas flow, patient's lung compliance and PIP settings control the inspiratory time and tidal volume. Then adjust rate dial to obtain desired breathing rate.

TABLE 1- ESTIMATED TIDAL VOLUME (ML) DELIVERED AT VARIOUS FLOW (LPM) AND INSPIRATORY TIME (SECONDS)

Flow (LPM)	Inspiratory Time (seconds)					
	0.5	1	1.5	2	2.5	3
15	125	250	375	500	625	750
20	167	333	500	667	833	1000
25	208	417	625	833	1042	1250
30	250	500	750	1000	1250	1500
35	292	583	875	1167	1458	1750
40	333	667	1000	1333	1667	2000

<sup>[1]</sup> Standard Specification for Minimum Performance and Safety Requirements for Resuscitators Intended for Use with Humans, ASTM Designation: F 920 - 93.

## VAR User's Guide

---

### II. Clinical Considerations (continued)

The VAR runs on a continuous gas flow (inspiratory flow) of 15 to 40 L/min depending on patients' inspiratory flow demand. When connected to a 50 PSIG gas source, the VAR will automatically deliver 40 L/min (667 mL/second) per ASTM guideline <sup>[1]</sup>. Delivered tidal volume may be determined by multiplying the flow in mL/second and the inspiratory time in second, or by using the estimated tidal volume table.

The rate dial controls exhalation time ( $t_{exhl}$ ), and when dialed down enough will cause the VAR to stop cycling automatically (infinite exhalation time). Under these circumstances the VAR is delivering pressure supported ventilatory support and the patient must trigger the VAR to begin subsequent full inhalations. If the patient is apneic or pressure control ventilation is desired, restart automatic cycling of the VAR by adjusting the rate dial counterclockwise until cycling begins again. Whenever the VAR stops cycling, the first step, in the absence of obvious clinical factors, is to check if it is in pressure support mode by rotating the rate dial counter clockwise (out). If rotating the rate dial counter clockwise substantially (3 or 4 turns) does not start automatic cycling, the patient's airway may be occluded or a very large leak exists.

Depending on the models for adult or pediatric patients, PIP may be adjusted from 10 and 50 cm H<sub>2</sub>O. The PEEP is intrinsic to the device which ranges from 2 to 9 centimeters and is directly proportional to the set PIP. Inspiratory time and rate are adjustable over a wide range. Changes in the PIP setting or flow will also affect the respiratory rate. It is important to check all settings when making a change to any of these three variables (flow, PIP and rate). For example: reducing the PIP setting may cause the VAR to go into spontaneous breathing mode. Adjust the rate dial out (counter-clockwise) to restart automatic cycling.

The VAR is equipped with an air entrainment valve which allows the patient to entrain additional air and respond to the demands of the patient. Patient entrainment of outside air is normally audibly detectable and the percent oxygen delivered to the patient will be reduced. Specific concentrations of oxygen may be delivered to the patient with the use of an oxygen blender.

Although the design of the modulator is similar to that of a pop-off valve and is inherently safe, the VAR is also equipped with a redundant pop-off valve that relieves pressure at 60 cm H<sub>2</sub>O. When the pop-off valve is activated, the pop-off valve piston will be seen to open slightly and excess pressure released.

---

## VAR User's Guide

---

### II. Clinical Considerations (continued)

Although peak pressures are listed on the side of the pressure dial, they are only approximate. Clinicians using the VAR are still required to use good clinical judgment and monitor the patient appropriately. A manometer may be connected between the modulator and the patient connector tee.

The VAR is pressure cycled on PEEP as well as PIP. In the pressure control mode there is no prolonged stage where the flow of exhalation gas stops for a significant duration of time (in the pressure support mode, exhalation time is determined by the patient). This occurs because the exhalation time is set with the rate dial by varying the exhalation resistance so that the patient just finishes exhalation with the beginning of the subsequent inhalation. The volume of gas with which the patient's lungs are inflated with when reaching PEEP is the same as with any other means of obtaining PEEP. As with all ventilatory support modes, short exhalation times on patients with high airway resistance may lead to gas trapping which is not detectable in the patient's external airways. Upon occlusion of the patient's airways, the VAR will stop cycling or may sometimes cycle rapidly.

The VAR will work with any mask that provides a good seal with the patient. All clinicians should receive adequate training with a VAR with mask prior to use. In the presence of a small leak, the VAR will still cycle between PIP and PEEP. Noticeable changes in the presence of a leak are increased inspiratory times and decreased expiratory times. The VAR works very well with an endotracheal tube.

Inhalation may be initiated by briefly removing the mask from the patient or briefly disconnecting the modulator from the patient adapter tee. In either event, inhalation begins because pressure has dropped down to PEEP and the VAR is pressure cycled.

Upon contamination of the VAR with vomitus, it may be cleared by disconnecting the modulator from the patient connector tee (see enclosed instructions) and tapping out vomitus on a hard surface. Additionally, if needed, the rate dial may also be removed to facilitate removal of vomitus from modulator. This operation should take less than 20 seconds, and in a lab setting has consistently been shown to take approximately 11 seconds. Alternatively, upon contamination with vomitus, the clinician may choose to discard the device and use a new one.


Inhalation and exhalation are audibly detectable and easily recognizable during operation of the VAR.

## VAR User's Guide

### II. Clinical Considerations (continued)

The VAR may be controlled remotely by connecting any length of 22-mm corrugated tubing between the patient connector tee and the modulator. The attached tubing will not increase the dead space, the modulator is an exhalation control valve, inspiratory gas is delivered through the patient connector tee.

The primary advantage of the VAR, as compared to manual resuscitators is the ability to deliver consistent, reliable, hands free resuscitation. Manual resuscitators may have adverse effects on patients as a result of inconsistent ventilation (see Clinical Reference in Section 5).



- 1 PIP – Set by **PIP DIAL**, control INSPIRATORY TIME ( $t_{insp}$ )
- 2 PEEP – Approximately 1/10<sup>th</sup> of PIP setting
- 3 INSPIRATORY FLOW RATE (Q) – Maximum 40 L/min (= 667 mL/sec)
- 4 INSPIRATORY TIME ( $t_{insp}$ ) – Time required to reach PIP
- 5 EXHALATION TIME ( $t_{exhl}$ ) – Time required to drop from PIP to PEEP
- 6 Tidal Volume =  $Q \times t_{insp}$
- 7 RESPIRATORY RATE (RR) =  $60 / (t_{insp} + t_{exhl})$
- 8 RATE DIAL – Set exhalation resistance and change RR

## VAR User's Guide

---

### III. Protocol: Setup Instructions VORTRAN Automatic Resuscitator

Policy Number:	Institution:	Department:
Date Adopted:	Dates Revised:	Dates Reviewed:
Approved by:	Name:	Title:

#### 1.0 POLICY STATEMENT:

This policy/protocol is intended for use with patients requiring short-term ventilatory support while being monitored by a clinician trained in the use of mechanical ventilation.

#### 2.0 PURPOSE:

To provide clinically appropriate recommendations and guidelines for the use of the VAR device, including clinical indications, device set-up, bedside application, potential hazards, and documentation.

#### 3.0 DESCRIPTION:

The VAR provides constant flow, pressure cycled ventilatory support in either pressure control or pressure support modes for breathing or non-breathing patients. The device includes a pulmonary modulator (an exhalation valve that opens at PIP and closes at PEEP) and a patient connector tee to supply gas flow, entrain additional air, and provide a redundant pressure pop-off valve for safety.

#### 4.0 PROCEDURE:

##### 4.1 INDICATIONS

- 4.1.1 Patients in need of emergency, short term, constant flow, pressure cycled ventilatory support.
- 4.1.2 Patients unable to maintain an adequate acid-base status during unassisted ventilation.

##### 4.2 CONTRAINDICATIONS

- 4.2.1 Patients requiring greater than 50 cm-H<sub>2</sub>O PIP.

##### 4.3 HAZARDS/PRECAUTIONS

- 4.3.1 The VAR should be used only by individuals who have adequate training in CPR techniques and the operation of gas-powered resuscitators.
- 4.3.2 Do not use grease or oil on the VAR for any reason.
- 4.3.3 Do not use the VAR in oxygen deficient atmospheres or near open flames.
- 4.3.4 Do not smoke while using the VAR or any other oxygen equipment.
- 4.3.5 Do not dismantle or attempt to remove any components other than those required for routine operations. Any tampering with the VAR may cause the unit to malfunction, and will automatically void the warranty.

## VAR User's Guide

---


### III. Protocol: VAR (continued)

#### 4.4 SET-UP INSTRUCTIONS

4.4.1 The VAR Models RT and RC are suitable for patients weighing over 88 pounds or 40 kilograms. For patients weighing 22 to 88 pounds, use our Model VAR Plus.

#### 4.4.2 **STEP [1]: SET FLOW**

Remove VAR from package and connect the supply tubing to either an appropriate cylinder or wall source. The VAR is designed to automatically deliver 40 liters per minute when connected directly to a 50 p.s.i.g. gas source.




**Note:** For better flow control, a flow meter capable of 40 liters per minute is preferred. The flow controls the inspiratory time – the higher the flow, the shorter the i-time; the lower the flow, the longer the i-time.

**Note:** If using an orifice-type flow regulator that is common to most cylinders, you will only be able to provide as much flow as the regulator indicates. If the regulator being used has a high flow port connection and you connect the VAR to this port, you will automatically get 40 liters per minute.


**Note:** The duration of an "E" cylinder when using a VAR will depend on the flow. An "E" cylinder contains 625 L of gas. At 40 L/min, 625 L will last up to 15 minutes; at 20 L/min, 625 L will last up to 30 minutes. 15 L/min orifice type flowmeters used on many "E" cylinders will not be able to deliver more than 15 L/min. When clinicians decide that 15 L/min is not sufficient flow, the VAR can be attached to a regulator that has a high flow port (50PSIG) to deliver 40 L/min. The length of use for various sizes of compressed oxygen tank (D, E, M & H) is a function of supplied oxygen flow from 6 to 40 L/min to VAR (see TABLE 1 - LENGTH OF USE FOR COMPRESSED OXYGEN TANKS).

**Note:** The VAR is completely gas driven, requiring no electrical power. The Model RT will deliver 100% oxygen to a patient when the VAR is supplied with 100% oxygen. If FiO<sub>2</sub> adjustment is desired, use oxygen blender at gas source.



**VAR User's Guide**

Note: The VAR Model RC and VAR-Plus Model PCM has a gas entrainment option for delivering FiO<sub>2</sub> at either 50 or 100%. If 100% is desired, connect one end of the supply tubing to the gas source and the other end to the green colored DISS connection marked "100%". The operating flow range is adjustable from 15 to 40 Liters per minute.



Note: If 50% FiO<sub>2</sub> is desired, remove the green "100%" adapter and connect the supply tubing to the gray colored entrainment adaptor on the patient tee. Set the oxygen supply flow from 6 to 15 Liters per minute. By entraining room air, the patient will receive flows of up to 40 liters per minute, significantly reducing the rate of oxygen consumption (refer to TABLE 2 - ENTRAINED FLOW TABLE).




TABLE 1 - LENGTH OF USE FOR COMPRESSED OXYGEN TANKS

Tank (Liters)	D	E	M	H
	387	622	3028	6905
Flow (LPM)	Length of use (minutes)			
6	65	100	500	1150
8	50	80	380	860
10	40	60	300	690
12	30	50	250	570
15	25	40	200	460
20	20	30	150	340
25	15	25	120	270
30	13	20	100	230
35	11	18	80	190
40	10	16	70	170

TABLE 2- - ENTRAINED FLOW TABLE

50% Connector		Total delivered flow (LPM)
Supply flow (LPM)	Entrained flow (LPM)	
6	14	20
8	17	25
10	20	30
12	23	35
15	25	40

## VAR User's Guide


---

### III. Protocol: VAR (continued)

#### 4.4 SET-UP INSTRUCTIONS

##### 4.4.3 **STEP [2]: SET PIP**

The second step in setting up the VAR is to set the patient's Peak Inspiratory Pressure, or PIP. Adjust the pressure dial to the desired setting.



- 4.4.4 Indicated peak pressure is printed on the PIP dial. Indicated pressures are approximate and depend on conditions and settings. Verify with a manometer by connecting between modulator and patient connector tee. I-time is counted off manually (1-1000, 2-1000, ) or with a watch.


Note: Positive End Expiratory Pressure, or PEEP, is intrinsic to the VAR. The PEEP ranges from 2 to 9 centimeters and is directly proportional to the set PIP.

Note: Indicated pressures are approximate and may vary depending on conditions and setting. Verify with a manometer

Note: Typical required supply pressure is 45 to 55 PSIG. Supply pressures from 39 to 80 PSIG may be used if the flow is adjusted to 40 L/min  $\pm$  10%. The VAR will deliver 40 L/min against a patient pressure of 20 to 40 cm-H<sub>2</sub>O when connected directly to a 50 PSIG source. Lower flows are obtainable with flowmeter adjustment. Use minimum flow of 20 L/min for best results.

##### 4.4.5 **STEP [3]: FUNCTION CHECK - CONNECT PATIENT**

Once your flow and pressure have been set, perform a function check on the unit before connecting it to the patient. This is accomplished by occluding the patient connection port and verifying that the modulator opens and the pressure does not exceed 60 cm-H<sub>2</sub>O.



- 4.4.6 **CONNECT PATIENT** - for use with a mask, clear mouth and airway of visible foreign bodies and use accepted techniques to ensure correct position of airways. Hold mask firmly against face while keeping head positioned. For use with endotracheal tube, connect endotracheal tube directly to patient adapter.

Note: It is very important to be trained in the correct application of the face mask before any attempt is made to use the VAR.


## VAR User's Guide

---

### III. Protocol: VAR (continued)

#### 4.4.7 **STEP [4] ADJUST RATE:**

Adjust Rate dial to achieve desired respiratory rate. This may only be accomplished after the patient is connected to the VAR.



Note: Remember, to slow down the patient's breathing rate, turn the rate dial clockwise. For a faster rate, turn the rate dial counter-clockwise.

Note: For a non-breathing patient, set the VAR in the Pressure Control Mode, which provides automatic cycling. Adjust the rate dial to the desired rate and count off the breathing rate manually.

If the VAR does not cycle, back out the rate dial until the device starts cycling again.

Note: If the patient is breathing spontaneously but needs ventilatory support, put the VAR in the Pressure Support Mode in which the patient triggers a breath. This is accomplished by adjusting the rate dial clockwise just until the VAR stops cycling automatically. Be sure to monitor the patient's effort with a manometer.

#### 4.4.8 **STEP [5] RE-ADJUST FLOW, PIP AND RATE:**

Observe the rise and fall of the chest corresponding to inhalation and exhalation of patient. Listen for expiratory flow from modulator. Listen to breath sounds of patient. There is no substitute for a good clinical assessment.

4.4.9 If the patient vomits, disconnect patient adapter from modulator and remove the rate dial if necessary. Tap out vomitus on a hard surface to dislodge and reassemble. Clear the patient's airway and reconnect. Clearing procedure should take less than 20 seconds. Check that inhalation and exhalation occur without obstruction.

4.4.10 The VAR is pressure limited and is equipped with a redundant pressure pop-off valve which will activate at 60 cm-H<sub>2</sub>O.

4.4.11 Changes in the patient's lung compliance will result in respiratory rate changes. In such an event, make the appropriate clinical changes.

4.4.12 If the patient draws air through the patient entrainment port, the oxygen concentration delivered to the patient may differ from the concentration at the gas inlet of the patient connector.

## VAR User's Guide

### IV. CAUTIONS AND WARNINGS

#### CAUTIONS

Federal law restricts the use of this device to sale by or on order of a physician (or properly licensed practitioner).

#### WARNINGS

1. The VORTRAN Automatic Resuscitator, VAR / VAR-Plus™ should be used only by individuals who have adequate training in CPR techniques and in the operation of gas powered resuscitators.
2. Do not use grease or oil on the VAR / VAR-Plus™ for any reason.
3. Spontaneously breathing patients may entrain ambient air.
4. Supply pressure of 39 to 80 PSIG must be adjustable to 40 L/min.
5. Redundant pop-off valve is set at 60 cm- H<sub>2</sub>O.
6. Do not use the VAR / VAR-Plus™ in oxygen deficient atmospheres or near open flames.
7. Do not smoke while using the VAR / VAR-Plus™ or any other oxygen equipment.
8. Do not dismantle or attempt to remove any components other than those required for routine operations. Any tampering with the VAR / VAR-Plus™ may cause the unit to malfunction and will automatically void the warranty.

#### PRECAUTION

1. Patients connected to this device are to be monitored continuously by persons having adequate training. Do not leave patients unattended.
2. When ventilating an intubated patient, higher pressure release settings may be required. Select a pressure setting of 35 cm H<sub>2</sub>O to start and adjust if necessary.
3. An audible, rapid clicking sound and rapid movement of the piston in the modulator indicates airway obstruction. Clear the airway and resume the ventilation procedure.
4. Positive End Expiratory Pressure (PEEP) is intrinsic to this device. For VAR™, PEEP is usually 1/10th PIP and will range from 2 to 5 cm H<sub>2</sub>O; and VAR-Plus™ PEEP is usually 1/5th PIP and will range from 2 to 9 cm-H<sub>2</sub>O depending on pressure settings. Verify actual PEEP with a manometer.
5. For a minute ventilation of 10 L/min and an I:E ratio of 1:1, [a] at 100% FiO<sub>2</sub> setting - the VAR / VAR-Plus™ will operate for 30 minutes (± 10%) with an output and supply flowrate set at 20 L/min on an "E" cylinder volume of 625 liters. [b] at 50% FiO<sub>2</sub> setting - the VAR / VAR-Plus™ will operate for 80 minutes (± 10%) with an output flow rate of 20 L/min and supply flowrate set at 8 L/min on an "E" cylinder volume of 625 liters.
6. Please review and follow the instructions and observe the warnings before using the VAR / VAR-Plus™.
7. If the use or operation of the VAR / VAR-Plus™ is unclear, contact your distributor or dealer for clarification.
8. The VAR / VAR-Plus™ is a resuscitation management system and should not be used as an unattended automatic ventilator.

**VAR User's Guide**

---

**V. VAR COMPETENCY**

How to set up your ventilator dependent patient using the **VORTRAN Automatic Resuscitator (VAR) Model RC**, a fully automatic disposable ventilator that operates with compressed gas.

Objectives

1. To be able to set up the VAR.
  - a. Setting the required flow for FiO<sub>2</sub> 100% or 50%
  - b. Getting the Peak Inspiratory Pressure and PEEP from the manometer
  - c. Adjusting the respiratory rate
  - d. For non-breathing and spontaneous breathing patient
  
2. To be able to troubleshoot and correct any problem that may arise with the use of the VAR.
  - a. Gas consumption during use
  - b. When adjusting the rate dial on the VAR, it sometimes stops cycling. What is happening?

Troubleshooting

1. I can set a constant respiratory rate and tidal volume with the VAR.  
 True                       False
  
2. With the VAR, compliance has a direct effect on the respiratory rate and volumes being delivered to your patient.  
 True                       False

After completion of the VAR competency, the practitioner should be able to set up the VAR and troubleshoot problems that may arise with its use.

Name	Institution:
Department	Dates Completed:

## VAR User's Guide

---

### VI. FAQ (Frequently Asked Questions)

Question	Answer
1. Nomenclatures	E-time: Exhalation time in seconds I-time: Inspiratory time in seconds L/min: Flow of gas in liters per minutes Manometer: Pressure gauge PIP: Peak inspiratory pressure PEEP: Positive end-expiratory pressure VAR: VORTRAN Automatic Resuscitator
2. How does <b>VAR</b> function during inhalation and exhalation?	<p>The <b>VAR</b> is a small automatic gas-powered resuscitator intended to provide pressure-limited, flow controlled ventilatory support for short-term emergency ventilatory support for both breathing and non-breathing patients while being monitored by a clinician or trained operator. The VAR is a single patient, multiple use device.</p> <ul style="list-style-type: none"> <li>• During inhalation, exhalation will not start until the desired peak inspiratory pressure (PIP) is reached.</li> <li>• During exhalation, inhalation will not begin until pressure drops to the controlled positive end-expiratory pressure (PEEP).</li> </ul>
3. What is the definition of "pressure control" mode (mandatory breathing) and "pressure support" mode (assisted breathing) when used with <b>VAR</b> ?	<p>For non-breathing patients - the mode of ventilation is called pressure control because no effort is required by the patient to initiate inhalation (mandatory breathing).</p> <p>For patients taking spontaneous breaths requiring assisted breathing - If the rate dial has been adjusted to a position that the continuous flow of gas does create more pressure than the set PEEP, then the <b>VAR</b> will not go into inhalation until the patient draws the baseline pressure down to PEEP. Again, because the <b>VAR</b> is cycled on both set PIP and PEEP, inhalation will not start until pressure reaches the set PEEP value. This mode of ventilation is called pressure support because the <b>VAR</b> only delivers ventilatory support when initiated by the patient.</p>
4. How do I set the <b>VAR</b> in pressure control or pressure support mode?	<p>Which mode the <b>VAR</b> is in is simply a function of where the rate dial has been adjusted. Turn the rate dial clockwise until in pressure support (assisted breathing), for pressure control (mandatory breathing) turn the rate dial counter clockwise.</p>
5. What does rate dial do?	<p>The rate dial is a variable resistor which controls the rate at which gas may escape. When the rate dial has been adjusted to a position that the continuous flow of gas does not create more pressure than the set PEEP (set PEEP is 1/10<sup>th</sup> of the set PIP, (for example: PIP=40 cm-H<sub>2</sub>O, PEEP=5 cm-H<sub>2</sub>O) then upon completion of exhalation the <b>VAR</b> will automatically cycle into inhalation because the <b>VAR</b> cycles on both the set PIP and PEEP.</p>
6. Does the gas supplied to the VAR flow continuously during exhalation and inhalation?	<p>Yes</p>

## VAR User's Guide

---

### VI. FAQ (Frequently Asked Questions)

Question	Answer
7. Does the <b>VAR</b> work with a mask or an endotracheal tube (Combitube® Dual Lumen Airway)?	The <b>VAR</b> works well with most endotracheal tubes or masks. When the <b>VAR</b> is used with a mask, clinicians must have proper training to be aware of the increased mortality associated with aspiration when vomitus occurs. If there is a small leak around the mask, the <b>VAR</b> will still cycle between PIP and PEEP, but inspiratory times will increase and expiratory times will decrease. In the event of a larger leak, the <b>VAR</b> will stop cycling because it cannot compensate for the leak.
8. What is the sensitivity or pressure drop required to trigger the <b>VAR</b> into inhalation?	The <b>VAR</b> is pressure cycled on PIP and PEEP. Therefore, as soon as the patient's pressure drops to PEEP, inhalation will start whether this occurs because exhalation has been completed or the patient draws a breath. Compared to time cycled ventilators, the sensitivity would be zero in the pressure control mode. In the pressure support mode the sensitivity may be set as light as 1 cm H <sub>2</sub> O or less; therefore, the patient's work of breathing will be minimal. If greater effort by the patient is desired, then it may be increased by turning the rate dial clockwise. Be sure to use a manometer when performing this procedure.
9. When adjusting the rate dial on the <b>VAR</b> , it sometimes stops cycling. What is happening?	The <b>VAR</b> 's rate dial controls rate by controlling the exhalation time. Once the PIP and inspiratory flow (LPM) have been set, inspiratory time is also set. The only way to control respiratory rate is by controlling the exhalation time. In the pressure control mode this is done with the rate dial which is actually a variable flow resistor. Depending on the patient and flow conditions used, it is possible to set the rate dial so that the continuous flow of gas always creates more pressure across the variable flow resistance than what the modulator is set to cycle at for PEEP, which means that the <b>VAR</b> is currently in the pressure support mode. In this condition the patient's airway pressure will remain slightly above set PEEP just as in a variable resistance PEEP valve and the <b>VAR</b> will not cycle. When pressure control is the mode of ventilation (which is required), the situation is easily corrected by dialing out the rate dial (counterclockwise) until the <b>VAR</b> starts cycling, thus reducing the variable resistance so that the patient's pressure is allowed to drop below PEEP and cycle the modulator automatically. In the pressure support mode it is the patient who initiates inhalation by drawing the baseline pressure down to the set PEEP value. Therefore, if the <b>VAR</b> is not cycling, chances are that the patient is not spontaneously breathing or the rate dial has been adjusted too far down, creating a baseline pressure which is too high above the set PEEP value for the patient to be able to initiate inhalation (the sensitivity is too high). In either event, turn the rate dial counterclockwise until the sensitivity is low enough for the patient to trigger inhalation, or the <b>VAR</b> will go into the pressure control mode.

## VAR User's Guide

---

### VI. FAQ (Frequently Asked Questions)

Question	Answer
10. If I connect to a 15 L/min flowmeter and dial it all the way up, what flow will I get through the <b>VAR</b> ?	Orifice type flowmeters like those which are commonly used on "E" cylinders will flow a maximum of what is indicated on the gauge. Timeter and other flowmeters using a floating ball as an indication of flow are capable of being adjusted to flows above what is indicated. If connecting to a Timeter flowmeter and adjusting the dial all the way open, the float will be slightly above the 15 L/min flow mark but 40 L/min will actually be flowing through the <b>VAR</b> . As long as the hospital gas supply and cylinder regulators are adjusted to 50 PSIG, which is the standard, the flow going through the <b>VAR</b> will never exceed 40 L/min.
11. What kinds of compressed gas source can I use with <b>VAR</b> ?	You can use any breathing gas from the hospital wall outlet or gas cylinder.
12. All I have are 15 L/min orifice type flowmeters with my "E" cylinders. 15 L/min of inspiratory flow is not enough flow for my patient. What can I do?	Some cylinder regulators equipped with an orifice type 15 L/min flowmeter are also equipped with a high flow (power take-off) port. If you connect the VAR to this port, you will automatically get 40 L/min. If 40 L/min is too much flow or you don't have a high flow port, you will need to use a different flowmeter.
13. How long will my "E" cylinder last with the <b>VAR</b> ?	It depends on the flowrate. There are 625 L in an "E" cylinder so at 40 L/min it will last approximately 15 minutes; at 20 L/min it will last 30 minutes.
14. What is the FiO <sub>2</sub> delivered to my patient?	When operated with pure oxygen, the VAR delivers 100% O <sub>2</sub> when there is no dilution. An optional dilution system (VAR model RC and VAR-Plus Model PCM) can be used to deliver an air-O <sub>2</sub> mixture of 50% FiO <sub>2</sub> and extend oxygen cylinder functional time.
15. May I connect any DISS connector to the patient tee connector threaded gas inlet fitting?	Yes
16. How can I measure tidal volume when using the VAR?	Tidal volume may be estimated by using the tidal volume chart included with the instructions. The VAR runs on a continuous fixed flow rate of gas (inspiratory flow) of up to 40 L/min (667 mL/second) when connected to a 50 PSIG gas source with associated flowmeter and control valve. Tidal volume is the inspiratory time multiplied by the flow rate (example: 1 second i-time X 667 mL/second = 667 mL tidal volume).
17. PIP ranges are indicated on the pressure dial, but what is the expected PEEP?	PEEP setting on <b>VAR</b> is automatically set at about 1/10 <sup>th</sup> of the selected PIP and it is good clinical practice to use a manometer to verify any pressure setting. PIP indications on the pressure dial are approximate only and ranges between 15 and 50 cm H <sub>2</sub> O, and PEEP ranges between 2 and 9 cm H <sub>2</sub> O respectively.

## VAR User's Guide

---

### VI. FAQ (Frequently Asked Questions)

Question	Answer
18. How do I connect a pressure manometer to the <b>VAR</b> ?	A manometer may be connected to the VAR by placing a 22 mm fitting between the modulator and patient connector tee (see enclosed instructions). Although the pressure dial indicates typical PIP and PEEP is 1/10th of PIP, a manometer is recommended because it provides valuable information to the clinician on what is occurring with the patient.
19. Is it possible to override the pop-off valve (high pressure relief valve)?	No, the <b>VAR</b> is a pressure cycled automatic resuscitator which has a maximum setting of 50 cm H <sub>2</sub> O. It includes an inspiratory pressure relief valve that opens automatically at approximately 60 cm H <sub>2</sub> O (preset and non-adjustable) and has a distinctive and easily recognized warning sound.
20. Can I deliver aerosol treatment while the patient is connected to the <b>VAR</b> ?	Yes. <b>NOTE:</b> Deposition of medicine residue may cause the <b>VAR</b> to stick if it dried for an extended period of time. Always perform a functional check per instructions before reconnecting the patient.
21. Is the <b>VAR</b> MRI compatible?	Yes. An extension kit with connecting tubing is available for MRI patient set up (refer to Technical Bulletin 082098 - MRI Compatibility).
22. Can I do CPR (closed-chest compression) with conventional automatic gas-powered resuscitators?	Yes. The cardiopulmonary resuscitation (CPR) guidelines and American Society for Testing and Materials caution against the use of automatic gas-powered resuscitators during CPR closed chest compression because the compression process may interfere with lung ventilation and airway resistance may prevent adequate ventilation.
23. Can I use the <b>VAR</b> with CPR?	Yes. The <b>VAR</b> is ideal for use in CPR. Our studies (by Otto Raabe, Ph.D. et. al) have shown that the <b>VAR</b> is safer than manual resuscitation using a BVM. The <b>VAR</b> should not cause baro-trauma, as the unit will not exceed the set peak inspiratory pressure and will automatically cycle at the end of each compression. In the case of manual bagging, medical personnel must be careful not to bag and compress the patient simultaneously in order to avoid high PIP. Manual bagging can cause pressures that can exceed 60 cm-H <sub>2</sub> O.
24. Is the <b>VAR</b> safe when used with CPR?	Because the <b>VAR</b> responds to thoracic pressure variations, it appears to provide the maximum ventilation possible during closed chest compression and responds with a full inhalation at high flow rate as soon as the compression ends. Because of its audible and visual indications of inhalation-exhalation cycling, elevated airway resistance or low tidal volume is readily observed by the rescuer.
25. How about the caution statement?	These results suggest that there is no contraindication associated with performing CPR closed chest compression while utilizing the <b>VAR</b> as a ventilatory resuscitator. Further, the results suggest that such use would be beneficial. A revision of CPR guidelines and ASTM 920-93 should be considered.

## VAR User's Guide

---

### VI. FAQ (Frequently Asked Questions)

Question	Answer
26. What are some of the other commonly used devices for providing patient ventilatory support?	The VAR is classified as a "Automatic Pressure-cycled, gas-powered resuscitator" per ASTM resuscitator guideline (F920-93). There are "Operator-powered resuscitator" such as Bag-Valve-Mask (BVM), "Manually-cycled, gas-powered resuscitator" such as Demand Valves; "Automatic-time cycled, gas-powered resuscitator" and "Volume-cycled, gas-powered resuscitator" such as the emergency transport ventilators.
27. Why should I use VAR when I am use to BVM?	"Operator-powered resuscitator" – Bag-Valve-Mask (BVMs) are the most commonly used devices for emergency short term ventilator support. They are typically disposable and are used extensively in the pre-hospital and inter-hospital markets. Manual resuscitators are labor intensive and are unable to deliver consistent ventilatory support. When used with a mask or endotracheal tube, they require the clinician to use both hands. They do not require being connected to a gas supply to provide ventilatory support but are almost always used in conjunction with compressed oxygen to increase the patient's FiO <sub>2</sub> . Although they appear easy to use, many studies have shown that they all deliver insufficient tidal volume and often deliver respiratory rates which are too high, resulting in significant adverse effects on the patients (refer to Section 5 - Clinical References). Nevertheless, many clinicians, when questioned about the use of manual resuscitators, feel certain that they personally deliver a consistent tidal volume of 750 mL per breath and that the ventilatory support they deliver is superior because of the feel they get through their hands when squeezing the bag.
28. What does auto-matic-cycled resuscitator do?	There are many automatic resuscitators that are gas or battery powered and non-disposable. All require some type of regular cleaning and are sold with some type of associated disposable products. Most are constant flow, time cycled devices with no high pressure relief valve which puts the patient in danger of a pneumothorax if there is an unexpected decrease in lung compliance. These devices have no monitoring or alarm features and have a minimum list price of \$1,000.
29. What is the least expensive automatic-cycled resuscitator cost?	The Oxylator EM-100 is a gas powered automatic resuscitator which provides constant flow, pressure cycled ventilatory support just like the VAR. Unlike the VAR, the Oxylator EM-100 is relatively heavy, non-disposable, and is not equipped with a pop-off valve. It listed at \$700.
30. What are the advantages of the transport ventilators?	Transport ventilators are equipped with sophisticated monitoring and alarm functions. They are usually able to provide several modes of ventilatory support and provide more versatile ventilation than the VAR. Because they are very complicated, significant training is needed. They are used in conjunction with disposable products and can cost as much as \$2,000 to \$5,000.

**VAR User's Guide**

**VII. Clinical Reference**

**Dave Swift, RRT, RRCP - Senior Therapist Ottawa Hospital, Ottawa, Ontario, Canada, Preparing for Mass Casualties & Mechanical Ventilation Alternatives, presented at 48<sup>th</sup> AARC International Respiratory Congress in Tampa, FL, Oct 5-7, 2002.**

The March 1995 Tokyo, Japan terrorist attack using the nerve agent Sarin sounded a wakeup call to health care workers. The intentional release of this neurotoxin resulted in 11 dead and five thousand exhibiting toxic symptoms. The health care system was rapidly overwhelmed.<sup>1</sup>

The National Capital Region of Ottawa is home to embassies of many nations and is viewed as a very high risk for a terrorist attack. As the sole Respiratory Therapist representative on the Chemical, Biological, Radiation and Nuclear Committee, it became rapidly apparent that there was a serious discrepancy between the number of ventilators available and the actual ventilator resources available. This finite limit was determined to be both unacceptable and avoidable. To avoid compromising patient care a cost effective method for treating the largest number of patients had to be determined.

It was determined that a pneumatic, automatic resuscitator offered the best clinical options. As it was not dependent on a/c power, was highly portable and relatively easy to use it seemed the most appropriate, cost effective choice.

The units were tested using the following clinical simulations: increased resistance, decreased compliance, increased compliance and with an air leak present. All units performed as advertised when faced with increased compliance, with delivered volumes decreasing and rates increasing with increased resistance and compliance. Serious clinical problems would be encountered with air leaks present and would need prompt medical intervention. Although all three units performed as advertised, each unit had individual characteristics that would have to be evaluated by the potential user as suitable for their own clinical applications.

The Vortran Automatic Resuscitator offered the capabilities of managing the largest number of patients at the most financially responsible cost. In addition, the unit has the advantage of ease of use and that the equipment offered a simple solution to the handling of contaminated units from a biological or terrorism incident, it is disposable. The costs of the other units prohibited one time use and would result in a lengthy and expensive decontamination process, which might also pose a hazard to hospital staff charged with decontamination.

Characteristics:	<b>AMBUMATIC<sup>2</sup></b>	<b>GENESIS II<sup>3</sup></b>	<b>VORTAN AUTOMATIC RESUSCITATOR<sup>4</sup></b>
Patient Type	Pediatric (>3 yrs) & adult	Pediatric (>3 yrs) & adult	Pediatric (>3 yrs) & adult
Power source	pneumatic	pneumatic	pneumatic
Portability	<1.5 lb.	<2 lb.	<11lb
Pressure cycled	yes	no	yes
Volume cycled	yes	yes	no
Rates	12 or 20	8 -12	0->40
Antisuffocation valve	yes	yes	yes
Pressure relief	yes	yes	yes
Pressure monitoring	yes (optional)	no	yes (optional)
Alarms	audible blowoff	audible blowoff	audible blowoff
FiO2 control	60 or 100%	100%	50 or 100%
PEEP	intrinsic	intrinsic	intrinsic
Single/multiple use	multiple pts	multiple pts	single
Cost CDN(0.62US\$)	>\$500	<\$400	<\$45
Replacement parts required & CT scan/ MRI compatibility	yes (valves, etc) Not certified for CT Scan or MRI use	no Not certified for CT Scan or MRI use	no Certified for CT Scan or MRI use

Characteristics Required In A Mass Casualty Ventilator/Resucitator:

<sup>1</sup> Brackett D.W., Holy Terror, Armageddon in Tokyo, New York:Weatherhill, Inc. 1996

<sup>2</sup> Ambu matic , Manufacturer: Ambu Inc. Linthicum, MD, USA

<sup>3</sup> GenesisII, Manufacturer :O2 Systems Inc., Mississauga, Ontario, Canada

<sup>4</sup> VAR (Resp. Tech Pro), Manufacturer: Vortran Medical Technology, Sacramento, California, USA

## VAR User's Guide

### VII. Clinical Reference

Steven J. Weiss, Todd Filbrun, Chad Augustin, Ray Jones and Amy Ernst . UC Davis Medical Center: Sacramento, CA, Sacramento City Fire/EMS: Sacramento, CA. ABSTRACT: An Automatic Transport Ventilator (ATV) vs. Bag Valve Mask (BVM) for Ventilation during EMS Transport. Academic Emergency Medicine Volume 11, Number 5 592, May 2004 .

**Abstract:**

**OBJECTIVES:** The hypothesis of this study was that paramedics (EMTPs) perceived that use of automatic transport ventilator (ATV) was better than BVM for managing ventilation during patient transport.

**METHODS:** ATVs and BVMs were randomized on 5 City Fire Department Paramedic Units. At the conclusion of each patient transport, using a 5-point Likert scale, EMTPs rated the modality used (ATV vs. BVM) on ease of use, time of setup, expedition of transport, additional tasks completed, documentation, overall patient care, and patient comfort. Pulse, oxygen saturation, respiratory rate, and end tidal CO<sub>2</sub> were collected every 5 seconds. Statistical analysis was performed on results of the Likert scale using a Mann-Whitney U rank sum test. Results were significant if  $p < 0.05$ . The power of the study was 80 percent to show a difference of 1.0 on the Likert scale.

**RESULTS:** 28 patients were entered into the study, 14 BVM and 14 ATV. The reason for device use was assisted ventilation in 7/28 (25%) cases and CPR in 21/28 (75%) cases. There were no significant differences in the EMS perceptions of ease of use ( $p = 0.08$ ), time of setup ( $p = 0.14$ ), expedition of transport ( $p = 0.27$ ), or overall patient care ( $p = 0.59$ ). There were significant differences in favor of the ATV in ability to accomplish additional tasks ( $p = 0.01$ ), ability to document ( $p = 0.04$ ), and ability to provide patient care ( $p = 0.03$ ). The data collector stored ongoing physiologic data on 15/28 (54%) patients during EMS transport.

**CONCLUSIONS:** EMTPs perceived that they were able to accomplish more tasks, document more completely, and provide better patient care with the use of the ATV. The data collector time marked data and stored the data for subsequent retrieval in a majority of cases.

*Nates, Joseph L. MD, FCCM* **Combined external and internal hospital disaster: Impact and response in a Houston trauma center intensive care unit** \*. Critical Care Medicine. 32(3):686-690, March 2004.

**Abstract:**

**Objective:** To increase awareness of specific risks to healthcare systems during a natural or civil disaster. We describe the catastrophic disruption of essential services and the point-by-point response to the crisis in a major medical center.

**Design:** Case report, review of the literature, and discussion.

**Setting:** A 28-bed intensive care unit in a level I trauma center in the largest medical center in the world.

**Case:** In June 2001, tropical storm Allison caused >3 feet of rainfall and catastrophic flooding in Houston, TX. Memorial Hermann Hospital, one of only two level I trauma centers in the community, lost electrical power, communications systems, running water, and internal transportation. All essential hospital services were rendered nonfunctional. Life-saving equipment such as ventilators, infusion pumps, and monitors became useless. Patients were triaged to other medical facilities based on acuity using ground and air ambulances. No patients died as result of the internal disaster.

**Conclusion:** Adequate training, teamwork, communication, coordination with other healthcare professionals, and strong leadership are essential during a crisis. Electricity is vital when delivering care in today's healthcare system, which depends on advanced technology. It is imperative that hospitals take the necessary measures to preserve electrical power at all times. Hospitals should have battery-operated internal and external communication systems readily available in the event of a widespread disaster and communication outage. Critical services such as pharmacy, laboratories, blood bank, and central supply rooms should be located at sites more secure than the ground floors, and these services should be prepared for more extensive performances. Contingency plans to maintain protected water supplies and available emergency kits with batteries, flashlights, two-way radios, and a nonelectronic emergency system for patient identification are also very important. Rapid adaptation to unexpected adverse conditions is critical to the successful implementation of any disaster plan.

## VAR User's Guide

### VII. Clinical Reference

Otto G. Raabe, Ph.D. and Mario Romano, RCP, Comparison of RespirTech PRO<sup>®</sup> and Ambu<sup>®</sup> SPUR Resuscitators During Simulated CPR.

**BACKGROUND:** The cardiopulmonary resuscitation (CPR) guidelines and American Society for Testing and Materials warn against the use of automatic pulmonary resuscitators during CPR closed chest compression because the compression process may interfere with lung ventilation and airway resistance may prevent adequate ventilation. However, appropriately designed pressure-cycled, pressure-controlled (rather than pressure-cycled, time-controlled) mechanical ventilators should be able to automatically respond to pulmonary pressure changes to provide air or oxygen to the lung at high flow rate upon demand and alert the rescuer of ventilatory problems. This evaluation was conducted to investigate ventilatory factors associated with the use of either the portable RespirTech PRO<sup>™</sup> (RTP) gas-powered automatic resuscitator or a typical manually operated self-inflating bag-valve resuscitator.

**METHODS:** Thirty tests, 17 with the RTP and 13 with the bag resuscitator, were conducted using the resuscitator connected to a commercial test lung modified for automatic simulated chest-compression following standard compression rates as timed with an electronic metronome. The test system was designed to be totally mechanical to avoid operator effects.

**RESULTS:** Both resuscitators provided appropriate ventilation without excessive lung pressures following the chosen 5:1 compression-ventilation ratio. Overall, the RTP (at 25 L/min) and bag-valve resuscitator minute ventilation values were about the same with means of  $6.3 \pm 0.5$  SE liters and  $6.2 \pm 0.6$  SE liters, respectively. The RTP automatically responded to pulmonary pressure variations, rapidly delivering short breaths between compressions and a full inhalation during the pause without serious pressure extremes. The highest observed intrapulmonary pressures ( $>80$  cm H<sub>2</sub>O) occurred with the bag-valve resuscitator operated during uninterrupted ("seamless") chest compressions without inhalation synchronization.

**DISCUSSION:** Both devices worked well following the standard protocol for CPR. Because the RTP inhalation-exhalation cycling is visually and audibly obvious, indications of possible airway resistance or low tidal volume are readily observed by the rescuer.

**CONCLUSIONS:** The RTP may be used safely as an automatic resuscitator during CPR. Revision of CPR guidelines and ASTM 920-93 for use of pressure-controlled resuscitators should be considered.

Michael Rossini, M.D., Barry Hickerson, EMT-P, Preliminary Evaluation of a Lightweight, Disposable Emergency Transport Ventilator in the Aeromedical Setting

**Introduction:** Recent evidence suggests patients receiving pre-hospital ventilation benefit from the use of emergency transport ventilators (ETV). This evidence is supported by the fact manual ventilation using a bag-valve-mask type device has substantial variations in rate and volume. These variations occur during initial treatment and transport even by well-trained EMS crews. Proper tidal volume, airway pressures and respiratory rate are critical components of emergency ventilatory support and variations can impact mortality and morbidity on a wide range of patients suffering from illness or injury.

**Methods:** The purpose of the evaluation was to determine the practicality and ease of use of a new ETV, the "RespirTech PRO" manufactured by Vortran Medical Technology 1, Inc. and identify any shortcomings during the initial phases of patient treatment, transport and emergency room care. The ETV was placed into service on our single BK 117-B2 hospital-based helicopter program. A Registered Nurse and Licensed Paramedic staff Air Med Team, which is based in Modesto, California. The majority of scene transports are flown to our base hospital, Doctors Medical Center also in Modesto, California. We gathered data on 12 patients from October 1999 to July 2000 that received ventilatory support from the ETV. Vitals signs during and post transport, arterial blood gases post transport and subjective data regarding ease of use, set-up and controls were gathered on all 12 patients.

**Results:** Twelve patients received on-scene and in-flight ventilatory support from the ETV without complications. All 12 adult patients were intubated by ground EMS personnel or the Air Med Team and placed on the ETV. The two manual setting, pressure and rate were set without difficulty and facilitated by the use of continuous end-tidal CO<sub>2</sub> monitoring. The oxygen source for the ventilator was a 15-25 liter per minute fitting that allowed operation without difficulty in all 12 cases. Blood gas analysis and review of vital signs during and post transport indicated all patients had been adequately ventilated during initial treatment and transport.

**Conclusion:** The RespirTech PRO proved to be an easy-to-use and reliable ETV that lends itself to a range of patients requiring prehospital ventilation. Ventilation is a key factor in the outcome of many types of injury and illness and this ETV should be considered for on-scene or transport use in a variety of prehospital settings.

## VAR User's Guide

### VII. Clinical Reference

Mario Romano, RCP, Otto G. Raabe, Ph.D., William Walby, MS and Timothy E. Albertson, MD, Ph.D., The Stability of Arterial Blood Gases During Transportation of Patients Using the RespirTech PRO™, American Journal of Emergency Medicine, May 2000.

**PURPOSE:** The transportation of critically ill patients requiring mechanical ventilation is recognized as a high risk and expensive procedure. Approaches have included using manual bag-type valve resuscitators and expensive portable transport ventilators. This study evaluated the effectiveness of the inexpensive portable RespirTech PRO (RTP) gas-powered automatic resuscitator during intra-hospital transport of critically ill mechanically ventilated patients.

**BASIC PROCEDURES:** Twenty medical intensive patients on stable mechanical ventilator settings had arterial blood gas and vital sign determination before routine transport out of the intensive care unit. Repeat measurements were made during transport approximately 30 minutes after being placed on the RTP portable pressure-cycled automatic resuscitator using an FiO<sub>2</sub> of 100%.

**MAIN FINDINGS:** During use of the RTP for transport, there were no statistically significant variations observed in mean arterial blood pressure [82 ± 11 SD (range 65-112) mm Hg before transport versus 85 ± 14 SD (range 59-110) mm Hg during transport], heart rate [94 ± 16 SD (range 74-127) beats/min before versus 96 ± 17 SD (range 69-132) beats/min during], arterial pH [7.41 ± 0.07 SD (range 7.31-7.58) before versus 7.42 ± 0.05 SD (range 7.37-7.52) during], and PaCO<sub>2</sub> [43 ± 10 SD (range 26-65) mm Hg versus during 43 ± 10 SD (range 27-61) mm Hg during]. Because the FiO<sub>2</sub> before transport was 63 ± 26 SD (range 30-100%) versus 100% during transport using the RTP, the mean PaO<sub>2</sub> was significantly increased from 124 ± 86 SD (range 57-367) mm Hg before transport to 297 ± 168 SD (range 65-537) mm Hg during (P<0.001). No transportation associated clinical adverse events were noted.

**DISCUSSION:** Several previous investigations have shown that portable ventilators are safe and effective in intra-hospital transport of mechanically ventilated patients. This study demonstrated that the portable pressure-cycled RTP also allows safe transportation of mechanically ventilated ICU patients. By analogy, the RTP is potentially useful as an automatic resuscitator for emergency medical care.

**PRINCIPAL CONCLUSION:** This RTP is a disposable resuscitator/ventilator device that provides an inexpensive alternative for transporting ventilator-dependent patients

Mario Romano, RCP, COMPATIBILITY OF THE RESPIRTECH PRO IN THE MRI UNIT, presented at 46th AARC International Respiratory Congress in Cincinnati, OH, October 7-10, 2000.

**BACKGROUND:** Because most medical facilities do not have MRI compatible ventilators, MRI studies on intubated patients are frequently delayed until the patient is extubated. Although there are mechanical ventilators that are MRI compatible, the cost for purchasing them for MRI use only is impractical, especially in light of the limited number of intubated patients needing an MRI. This paper examines the RespirTech PRO, a single patient use fully automatic resuscitator, and how it functioned during an MRI study in a General Electric 1.5 MRI unit.

**METHODS:** One clinically stable 72-year old male patient in need of an MRI of his head was placed on the automatic resuscitator with extension kit. The patient was set in a control mode of 16 BPM with a ventilating pressure of 25 cm-H<sub>2</sub>O and a liter flow of 40 LPM at a FiO<sub>2</sub> of 100%.

The patient was placed in a General Electric 1.5 MRI unit, and the device functioned without incident. No attraction to the magnet was noted. Image artifact was minimal and was limited to the patient tee area, allowing for a clear picture of the head. The patient tolerated ventilation well, and his vital signs are summarized in the graph below.

<b>RESULTS:</b> Patient Vital Signs:	Tx	HR	BP	O <sub>2</sub> Sat.	FiO <sub>2</sub> Set
Pre MRI	85	98 / 51	98%	100%	
During MRI	77	102 / 54	96%	100%	

**DISCUSSION:** No significant changes in vital signs or O<sub>2</sub> saturation were noted with the use of the automatic resuscitator. The patient appeared to tolerate the procedure with no adverse affects. No attraction to the MRI magnet was noted and artifact was limited to the patient tee area.

**Conclusions:** The RespirTech PRO can be a safe and cost effective ventilator for use in the MRI room without the need to purchase capital equipment. More experience with the use of this automatic resuscitator in transporting patients to other areas of the hospital can establish it as a safe and cost effective transport.

---

**VAR User's Guide**

---

**VIII. Coding information**

**HCPCS - HCFA (Health Care Financing Administration)  
Common Procedure Coding System**

PRODUCT ----- VORTRAN Automatic Resuscitator

CODE ----- K0533 / E0471

DESCRIPTION - Respiratory assist device, bi-level pressure capability,  
with backup rate feature, used with noninvasive inter-  
face, e.g. nasal or facial mask (intermittent assist device  
with continuous positive airway pressure device)

INSTRUCTIONS Coverage issue, CIM 60-9

**CPT - Current Procedure Terminology  
(American Medical Association)**

PRODUCT ----- VORTRAN Automatic Resuscitator

CODE ----- 94656

DESCRIPTION - Ventilation assist and management, initiation of pres-  
sure or volume preset ventilators for assisted or con-  
trolled breathing; first day.

CODE ----- 94657

DESCRIPTION - Subsequent

---

## VAR User's Guide









---

### IX. Troubleshooting

	<b>Cause / Action</b>
• VAR stops cycling	<ol style="list-style-type: none"> <li>1. Leak in circuit / look for pressure leak at gas connections, gas supply, patient airway, ET tube cuff, etc.</li> <li>2. Compliance change / change in patient's lung compliance may effect breathing rate, adjust rate dial accordingly.</li> <li>3. Increased gas flow / an increase in supply gas flow will elevate PEEP. Adjust rate dial out accordingly.</li> <li>4. Out of gas / replace oxygen cylinder or connect to other gas source.</li> <li>5. Assisted breathing mode / VAR is waiting for patient to trigger the inspiration.</li> </ol>
• Breathing (exhalation) rate is too fast	<ol style="list-style-type: none"> <li>1. High supply flow / reduce supply gas flow.</li> <li>2. Low PIP setting / increase PIP as needed.</li> <li>3. Compliance too high (stiff lung) / change mode of ventilation.</li> <li>4. Low exhalation resistance / set rate dial down (clockwise)</li> </ol>
• Breathing (exhalation) rate is too slow	<ol style="list-style-type: none"> <li>1. Low supply flow / increase supply gas flow.</li> <li>2. High PIP setting / decrease PIP as needed.</li> <li>3. Compliance too low (soft lung) / change mode of ventilation.</li> <li>4. High exhalation resistance / adjust rate dial out (counter-clockwise).</li> </ol>
• Inspiratory time is too long	<ol style="list-style-type: none"> <li>1. Supply gas flow too low / increase supply gas flow.</li> <li>2. High PIP setting / lower PIP setting as appropriate.</li> <li>3. Compliance too low / change mode of ventilation.</li> </ol>
• Inspiratory time is too short	<ol style="list-style-type: none"> <li>1. Supply gas flow too high / decrease supply gas flow.</li> <li>2. Low PIP setting / increase PIP setting as appropriate.</li> <li>3. Compliance too high / change mode of ventilation.</li> </ol>
• Reading on pressure manometer increase or decrease	<ol style="list-style-type: none"> <li>1. PIP increase / look for airway occlusion or kinked ET tube, change in patient's lung compliance, verify PIP setting.</li> <li>2. PIP decrease / look for change in patient's lung compliance, verify PIP setting on VAR.</li> </ol>
• Flow on my regulator does not go to 40 LPM	<ol style="list-style-type: none"> <li>1. VAR can operate at flows as low as 15 LPM</li> <li>2. Maximum flow is 40 LPM per ASTM guideline</li> <li>3. If you have a 15-16 LPM flowmeter and is connect to 50 PSIG piped in gas source, the VAR will flow limited at 40 LPM when you flush open the flowmeter.</li> </ol>
• Maximum flow 40 LPM is not enough	<ol style="list-style-type: none"> <li>1. Maximum flow is 40 LPM per ASTM guideline.</li> <li>2. Patient can entrain additional room air through the one-way valve to meet their inspiratory flow demand for spontaneously breathing patient.</li> </ol>
• What is the tidal volume delivered	<ol style="list-style-type: none"> <li>1. Tidal volume delivered is a function of flow rate (liters per minute) over inspiratory time (seconds).</li> <li>2. VAR delivers gas until the set PIP is reached.</li> </ol>






**VAR User's Guide**

**X. Ordering information**

Product Description	Order No	Case Quantity	Accessories					
			O <sub>2</sub> tubing	Manometer	Flex hose	22 mm tubing	Entrainment	
<b>Model VAR-Plus for patients weighing between 10 to 40 kilograms</b>								
<b>VAR-Plus Model PCM</b>	PCM- 5011	10	7'	√	√		√	
<b>VAR Models RT and RC for patients weighing over 40 kilograms</b>								
<b>VAR Model RT</b>	RT- 4000	10	7'		√			
VAR RT w/Manometer	RTM- 4001	10	7'	√	√			
VAR RT w/Extension	RTE- 4002	10	20'	√	√	√		
VAR RTM 4 pack	RTF- 4003	4	7'	√	√			
<b>VAR Model RC</b>	RC- 4001	10	7'		√		√	
VAR RC w/Manometer	RCM- 4011	10	7'	√	√		√	
VAR RC w/Extension	RCE- 4012	10	20'	√	√	√	√	
Regulator w/ 4 RCM	RCK- 4050	4	7'	√	√		√	
Regulator, Oxygen	RG- 2162	1						
Extension Kit only	EXK- 4040	10	20'			√		

**VAR User's Guide**

**XI. Quick Guide**

<p><b>1</b></p>	<p><b>Set flow</b></p> <p>Step 1, set the flow to approximately 25 liters per minute with a flow meter from a suitable gas source.</p>	
<p><b>2</b></p>	<p><b>Set PIP</b></p> <p>Step 2, set the PIP with the pressure dial. Each VAR is factory tested and pre-set to 25 centimeters of water pressure.</p>	
<p><b>3</b></p>	<p><b>Function check. Connect to patient</b></p> <p>Step 3, before connecting to the patient, perform a function check. Connect the VAR to the patient and then verify the PIP with a pressure manometer.</p>	
<p><b>4</b></p>	<p><b>Adjust breathing rate</b></p> <p>Step 4, adjust the breathing rate with the rate dial.</p>	
<p><b>5</b></p>	<p><b>Adjust Flow, PIP and Rate</b></p> <p>Step 5, re-adjust flow, PIP and rate for the patient's clinical situation. There is no substitute for a good clinical assessment.</p>	
<p><b>*</b></p>	<p><b>Use of extension tubing in MRI or CT</b></p> <p>NOTE: Position modulator away from patient tee using extension tubing and verify patient cycling whenever modulator was repositioned.</p>	

This Quick Guide is intended to help you gain a general understanding of the VAR product. Please be certain to read, understand, and follow the information listed in this User's Guide before using this product.

© 2005 - VORTRAN Medical Technology 1, Inc., Sacramento, California U.S.A.